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## Report HTA e HTA-like

1. Ontario Health. *Use of B-type natriuretic peptide (BNP) and N-terminal proBNP (NT-proBNP) as diagnostic tests in adults with suspected heart failure: a health technology assessment*. Ont Health Technol Assess Ser. 2021

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2. Hendricks S, Dykun I, Balcer B, Totzeck M, Rassaf T, Mahabadi AA., *Higher BNP/NT-pro BNP levels stratify prognosis equally well in patients with and without heart failure: a meta-analysis*, ESC Heart Fail 2022, doi: 10.1002/ehf2.14019
3. Rimmelzwaal S, van Ballegooijen AJ, Schoonmade LJ, Dal Canto E, Handoko ML, Henkens MTHM, van Empel V, Heymans SRB, Beulens JWJ., *Natriuretic peptides for the detection of diastolic dysfunction and heart failure with preserved ejection fraction-a systematic review and meta-analysis*, BMC Med 2020, doi: 10.1186/s12916-020-01764-x
4. Ibrahim M, Ahmad J, Abbas M, Zainullah, Umar Z, Nasir M, Zain K, Ahmad J, Arshad S, Bashir A, Ullah S, Ahmad Z, Safdar S., *The Role of N-terminal Pro-B-Type Natriuretic Peptide, Troponins, and D-dimer in Acute Cardio-Respiratory Syndromes: A Multi-specialty Systematic Review*, Cureus 2025, doi: 10.7759/cureus.84460
5. Vulciu PA, Pilat L, Mot MD, Barata PI, Mikos ID, Ioana MRS, Alexandru A, Ivan CS, Varga NI, Mladin NC, Puschita M., *Diagnostic and Prognostic Potential of Tetranectin in Heart Failure and Cardiovascular Disease: A Systematic Review*, Med Sci (Basel) 2025, doi: 10.3390/medsci13040206
6. van Dijk SHB, Brusse-Keizer MGJ, Bucsán CC, Ploumen EH, van Beurden WJC, van der Palen J, Doggen CJM, Lenferink A., *Lack of Evidence Regarding Markers Identifying Acute Heart Failure in Patients with COPD: An AI-Supported Systematic Review*, Int J Chron Obstruct Pulmon Dis 2024, doi: 10.2147/COPD.S437899
7. Rao RA, Varghese SS, Ansari F, Rao A, Meng E, El-Diasty M., *The Role of Natriuretic Peptides in Predicting Adverse Outcomes After Cardiac Surgery: An Updated Systematic Review*, Am J Cardiol 2024, doi: 10.1016/j.amjcard.2023.09.101
8. Janssen E, Jukema JW, Beeres SLMA, Schalijs MJ, Tops LF., *Prognostic Value of Natriuretic Peptides for All-Cause Mortality, Right Ventricular Failure, Major Adverse Events, and Myocardial Recovery in Advanced Heart Failure Patients Receiving a Left Ventricular Assist Device: A Systematic Review*, Front Cardiovasc Med 2021, doi: 10.3389/fcvm.2021.699492
9. Huixing L, Di F, Daoquan P., *Effect of Glucagon-like Peptide-1 Receptor Agonists on Prognosis of Heart Failure and Cardiac Function: A Systematic Review and Meta-analysis of Randomized Controlled Trials*, Clin Ther 2023, doi: 10.1016/j.clinthera.2022.12.006
10. Parvan R, Hosseinpour M, Moradi Y, Devaux Y, Cataliotti A, da Silva GJJ., *Diagnostic performance of microRNAs in the detection of heart failure with reduced or preserved ejection fraction: a systematic review and meta-analysis*, Eur J Heart Fail 2022, doi: 10.1002/ehf.2700
11. Yan H, Ma F, Zhang Y, Wang C, Qiu D, Zhou K, Hua Y, Li Y., *miRNAs as biomarkers for diagnosis of heart failure: A systematic review and meta-analysis*, Medicine (Baltimore) 2017, doi: 10.1097/MD.0000000000006825
12. van Doorn S, Geersing GJ, Kievit RF, van Mourik Y, Bertens LC, van Riet EES, Boonman-de Winter LJ, Moons KGM, Hoes AW, Rutten FH., *Opportunistic screening for heart failure with natriuretic peptides in patients with atrial fibrillation: a meta-analysis of individual participant data of four screening studies*, Heart 2018, doi: 10.1136/heartjnl-2017-312781

13. Henkens MTHM, Remmelzwaal S, Robinson EL, van Ballegooijen AJ, Barandiarán Aizpurua A, Verdonschot JAJ, Raafs AG, Weerts J, Hazebroek MR, Sanders-van Wijk S, Handoko ML, den Ruijter HM, Lam CSP, de Boer RA, Paulus WJ, van Empel VPM, Vos R, Brunner-La Rocca HP, Beulens JWJ, Heymans SRB., *Risk of bias in studies investigating novel diagnostic biomarkers for heart failure with preserved ejection fraction. A systematic review*, Eur J Heart Fail 2020, doi: 10.1002/ejhf.1944
14. Shi Y, Liu J, Liu C, Shuang X, Yang C, Qiao W, Dong G., *Diagnostic and prognostic value of serum soluble suppression of tumorigenicity-2 in heart failure with preserved ejection fraction: A systematic review and meta-analysis*, Front Cardiovasc Med 2022, doi: 10.3389/fcvm.2022.937291
15. Lee KK, Doudesis D, Anwar M, Astengo F, Chenevier-Gobeaux C, Claessens YE, Wussler D, Kozhuharov N, Strebel I, Sabti Z, deFilippi C, Seliger S, Moe G, Fernando C, Bayes-Genis A, van Kimmenade RRJ, Pinto Y, Gaggin HK, Wiemer JC, Möckel M, Rutten JHW, van den Meiracker AH, Gargani L, Pugliese NR, Pemberton C, Ibrahim I, Gegenhuber A, Mueller T, Neumaier M, Behnes M, Akin I, Bombelli M, Grassi G, Nazerian P, Albano G, Bahrman P, Newby DE, Japp AG, Tsanas A, Shah ASV, Richards AM, McMurray JJV, Mueller C, Januzzi JL, Mills NL; CoDE-HF investigators., *Development and validation of a decision support tool for the diagnosis of acute heart failure: systematic review, meta-analysis, and modelling study*, BMJ 2022, doi: 10.1136/bmj-2021-068424
16. Gallagher J, McCormack D, Zhou S, Ryan F, Watson C, McDonald K, Ledwidge MT., *A systematic review of clinical prediction rules for the diagnosis of chronic heart failure*, ESC Heart Fail 2019, doi: 10.1002/ehf2.12426
17. Ross L, Moxey J, Nikpour M., *Are troponin and B-type natriuretic peptides useful biomarkers for the diagnosis of systemic sclerosis heart involvement? A systematic literature review*, Semin Arthritis Rheum 2021, doi: 10.1016/j.semarthrit.2020.10.016
18. Ur Rehman M, Saeed H, Babar ZUD, Omer O, Tashfeen S., *Metabolomic Biomarkers in Heart Failure: A Systematic Review of Diagnostic and Prognostic Significance*, Cureus 2025, doi: 10.7759/cureus.97636
19. Goyder CR, Roalfe AK, Jones NR, Taylor KS, Plumptre CD, James O, Fanshawe TR, Hobbs FDR, Taylor CJ., *Diagnostic accuracy of natriuretic peptide screening for left ventricular systolic dysfunction in the community: systematic review and meta-analysis*, ESC Heart Fail 2023, doi: 10.1002/ehf2.14314
20. Homar V, Mirosevic S, Svab I, Lainscak M., *Natriuretic peptides for heart failure screening in nursing homes: a systematic review*, Heart Fail Rev 2021, doi: 10.1007/s10741-020-09944-w
21. Martindale JL, Wakai A, Collins SP, Levy PD, Diercks D, Hiestand BC, Fermann GJ, deSouza I, Sinert R., *Diagnosing Acute Heart Failure in the Emergency Department: A Systematic Review and Meta-analysis*, Acad Emerg Med 2016, doi: 10.1111/acem.12878
22. Ammar LA, Massoud GP, Chidiac C, Booz GW, Altara R, Zouein FA., *BNP and NT-proBNP as prognostic biomarkers for the prediction of adverse outcomes in HFpEF patients: A systematic review and meta-analysis*, Heart Fail Rev 2025, doi: 10.1007/s10741-024-10442-6
23. Singh G, Bamba H, Inban P, Chandrasekaran SH, Priyatha V, John J, Prajjwal P., *The role of biomarkers in the prognosis and risk stratification in heart failure: A systematic review*, Dis Mon 2024, doi: 10.1080/00015385.2024.2371628
24. El-Sherbini AH, Zucker J, Sartor C, Liblik K, Cherif IA, Fatehi Hassanabad A, Fernandez AL, El-Diasty M., *The Role of Natriuretic Peptides in Pericardial Fluid in Predicting Cardiovascular Disorders: A Systematic Review*, Cardiol Rev 2024, doi: 10.1097/CRD.0000000000000779
25. Chen H, Chhor M, Rayner BS, McGrath K, McClements L., *Evaluation of the diagnostic accuracy of current biomarkers in heart failure with preserved ejection fraction: A systematic review and meta-analysis*, Arch Cardiovasc Dis 2021, doi: 10.1016/j.acvd.2021.10.007
26. Taylor KS, Verbakel JY, Feakins BG, Price CP, Perera R, Bankhead C, Plüddemann A., *Diagnostic accuracy of point-of-care natriuretic peptide testing for chronic heart failure in ambulatory care: systematic review and meta-analysis*, BMJ 2018, doi: 10.1136/bmj.k1450
27. Dogheim GM, Amralla MT, Werida RH., *The clinical significance of neopterin and NT-pro BNP in chronic heart failure: a systematic review*, Acta Cardiol 2024, doi: 10.1080/00015385.2024.2371628
28. Arshad HS, Laraib F, Nikeeta F, Vandina F, Alias Gurya A, Farooq U, Ali S, Nisar L, Zamir MH, Akbar F., *Prognostic Utility of N-terminal Pro-B-type Natriuretic Peptide (NT-proBNP) for Predicting Hospital*

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29. Jacobson TA, Rahbari KJ, Schwartz WA, Bae Y, Zhang R, Nunes DA, Huang C, Issa RP, Smilowitz K, Yan LD, Hirschhorn LR, Khan SS, Huffman MD, Miller GE, Feinglass JM, McDade TW, Funk WE., *Dried Blood Spots to Assess Cardiovascular-Kidney-Metabolic Health*, J Am Heart Assoc 2025, doi: 10.1161/JAHA.124.037454
30. Aimo A, Januzzi JL Jr, Vergaro G, Ripoli A, Latini R, Masson S, Magnoli M, Anand IS, Cohn JN, Tavazzi L, Tognoni G, Gravning J, Ueland T, Nymo SH, Rocca HB, Bayes-Genis A, Lupón J, de Boer RA, Yoshihisa A, Takeishi Y, Egstrup M, Gustafsson I, Gaggin HK, Eggers KM, Huber K, Tentzeris I, Wilson Tang WH, Grodin JL, Passino C, Emdin M., *High-sensitivity troponin T, NT-proBNP and glomerular filtration rate: A multimarker strategy for risk stratification in chronic heart failure*, Int J Cardiol 2019, doi: 10.1016/j.ijcard.2018.10.079
31. Parvan R, Becker V, Hosseinpour M, Moradi Y, Louch WE, Cataliotti A, Devaux Y, Frisk M, Silva GJJ; AtheroNET COST Action CA21153., *Prognostic and predictive microRNA panels for heart failure patients with reduced or preserved ejection fraction: a meta-analysis of Kaplan-Meier-based individual patient data*, BMC Med 2025, doi: 10.1186/s12916-025-04238-0
32. Marx N, Federici M, Schütt K, Müller-Wieland D, Ajjan RA, Antunes MJ, Christodorescu RM, Crawford C, Di Angelantonio E, Eliasson B, Espinola-Klein C, Fauchier L, Halle M, Herrington WG, Kautzky-Willer A, Lambrinou E, Lesiak M, Lettino M, McGuire DK, Mullens W, Rocca B, Sattar N and ESC Scientific Document Group, *2023 ESC Guidelines for the management of cardiovascular disease in patients with diabetes*, doi: 10.1093/eurheartj/ehad192
33. Vergani M, Cannistraci R, Perseghin G, Ciardullo S, *Use of natriuretic peptides in screening, diagnosis, and monitoring of heart failure*, The Journal of AMD 2024, doi: 10.36171/jamd 24.27.1.8
34. Pop-Busui R, Repetto E, Baron J, Schumacher D, Vaduganathan M, Pandey A, *Screening Natriuretic Peptide Levels Predicts Heart Failure and Death in Individuals With Type 1 and Type 2 Diabetes Without Known Heart Failure*, Diabetes Care 2025, doi: 10.2337/dc25-1260
35. Kristensen SL, Mogensen UM, Jhund PS, Rørth R, Anand IS, Carson PE, Desai AS, Pitt B, Pfeffer MA, Solomon SD, Zile MR, Køber L, McMurray JJV., *N-Terminal Pro-B-Type Natriuretic Peptide Levels for Risk Prediction in Patients With Heart Failure and Preserved Ejection Fraction According to Atrial Fibrillation Status*, Circ Heart Fail 2019, doi: 10.1161/CIRCHEARTFAILURE.118.005766
36. Mahmood A, Dhall E, Primus CP, Gallagher A, Zakeri R, Mohammed SF, Chahal AA, Ricci F, Aung N, Khanji MY., *Heart failure with preserved ejection fraction management: a systematic review of clinical practice guidelines and recommendations*, Eur Heart J Qual Care Clin Outcomes 2024, doi: 10.1093/ehjqcco/qcae053
37. Emdin M, Clerico A, Clemenza F, Galvani M, Latini R, Masson S, Mulè P, Panteghini P, Valle R, Zaninotto M, Ganau A, Mariotti R, Volpe M, Aspromonte N, Cacciatore G, Cappelletti P, L'Abbate A, Miglio F, Ottani F, Pagani F, Passino C, Plebani M, Sarzani R, Zucchelli G, *Raccomandazioni sull'Impiego Clinico dei Peptidi Natriuretici Cardiaci*, Ital Heart J Suppl 2005
38. American Diabetes Association Professional Practice Committee for Diabetes, *Cardiovascular Disease and Risk Management: Standards of Care in Diabetes—2026*, Diabetes Care 2026;49(Supplement\_1), doi: 10.2337/dc26-S010
39. Castiglione V, Aimo A, Vergaro G, Saccaro L, Passino C, Emdin M. Biomarkers for the diagnosis and management of heart failure. Heart Fail Rev. 2022 Mar;27(2):625-643. doi: 10.1007/s10741-021-10105-w. Epub 2021 Apr 14. PMID: 33852110; PMCID: PMC8898236.
40. Bruno G, Landi A, Barutta F, Ghezzo G, Baldin C, Spadafora L, Schimmenti A, Prinzi T, Cavallo Perin P, Gruden G. N-terminal pro-brain natriuretic peptide is a stronger predictor of cardiovascular mortality than C-reactive protein and albumin excretion rate in elderly patients with type 2 diabetes: the Casale Monferrato population-based study. Diabetes Care. 2013 Sep;36(9):2677-82. doi: 10.2337/dc13-0353. Epub 2013 Apr 5. PMID: 23564918; PMCID: PMC3747916.

41. Hillis GS, Welsh P, Chalmers J, Perkovic V, Chow CK, Li Q, Jun M, Neal B, Zoungas S, Poulter N, Mancia G, Williams B, Sattar N, Woodward M. The relative and combined ability of high-sensitivity cardiac troponin T and N-terminal pro-B-type natriuretic peptide to predict cardiovascular events and death in patients with type 2 diabetes. *Diabetes Care*. 2014;37(1):295-303. doi: 10.2337/dc13-1165. Epub 2013 Oct 2. Erratum in: *Diabetes Care*. 2015 Jul;38(7):1393. doi: 10.2337/dc15-er07. PMID: 24089534.
42. Bidadkosh A, Lambooy SPH, Heerspink HJ, Pena MJ, Henning RH, Buikema H, Deelman LE. Predictive Properties of Biomarkers GDF-15, NTproBNP, and hs-TnT for Morbidity and Mortality in Patients With Type 2 Diabetes With Nephropathy. *Diabetes Care*. 2017 Jun;40(6):784-792. doi: 10.2337/dc16-2175. Epub 2017 Mar 24. PMID: 28341782.
43. Jarolim P, White WB, Cannon CP, Gao Q, Morrow DA. Serial Measurement of Natriuretic Peptides and Cardiovascular Outcomes in Patients With Type 2 Diabetes in the EXAMINE Trial. *Diabetes Care*. 2018 Jul;41(7):1510-1515. doi: 10.2337/dc18-0109. Epub 2018 May 30. PMID: 29848776.
44. Huelsmann M, Neuhold S, Resl M, Strunk G, Brath H, Francesconi C, Adlbrecht C, Prager R, Luger A, Pacher R, Clodi M. PONTIAC (NT-proBNP selected prevention of cardiac events in a population of diabetic patients without a history of cardiac disease): a prospective randomized controlled trial. *J Am Coll Cardiol*. 2013 Oct 8;62(15):1365-72. doi: 10.1016/j.jacc.2013.05.069. Epub 2013 Jun 27. PMID: 23810874.
45. Neuhold S, Resl M, Huelsmann M, Strunk G, Adlbrecht C, Rath C, Prager R, Luger A, Clodi M, Pacher R. Repeat measurements of glycated haemoglobin A(1c) and N-terminal pro-B-type natriuretic peptide: divergent behaviour in diabetes mellitus. *Eur J Clin Invest*. 2011 Dec;41(12):1292-8. doi: 10.1111/j.1365-2362.2011.02539.x. Epub 2011 May 25. PMID: 21615391.
46. Scirica BM, Braunwald E, Raz I, Cavender MA, Morrow DA, Jarolim P, Udell JA, Mosenzon O, Im K, Umez-Eronini AA, Pollack PS, Hirshberg B, Frederich R, Lewis BS, McGuire DK, Davidson J, Steg PG, Bhatt DL; SAVOR-TIMI 53 Steering Committee and Investigators\*. Heart failure, saxagliptin, and diabetes mellitus: observations from the SAVOR-TIMI 53 randomized trial. *Circulation*. 2014 Oct 28;130(18):1579-88. doi: 10.1161/CIRCULATIONAHA.114.010389. Epub 2014 Sep 4. Erratum in: *Circulation*. 2015 Oct 13;132(15):e198. doi: 10.1161/CIR.0000000000000330. PMID: 25189213.
47. Hamano K, Nakadaira I, Suzuki J, Gonai M. N-terminal fragment of probrain natriuretic peptide is associated with diabetes microvascular complications in type 2 diabetes. *Vasc Health Risk Manag*. 2014 Oct 3;10:585-9. doi: 10.2147/VHRM.S67753. Erratum in: *Vasc Health Risk Manag*. 2014;10:649. PMID: 25328404; PMCID: PMC4199566.
48. Zelniker TA, Morrow DA, Mosenzon O, Goodrich EL, Jarolim P, Murphy SA, Bhatt DL, Leiter LA, McGuire DK, Wilding J, Bode C, Lewis BS, Gause-Nilsson I, Langkilde AM, Fredriksson M, Raz I, Sabatine MS, Wiviott SD. Relationship between baseline cardiac biomarkers and cardiovascular death or hospitalization for heart failure with and without sodium-glucose co-transporter 2 inhibitor therapy in DECLARE-TIMI 58. *Eur J Heart Fail*. 2021 Jun;23(6):1026-1036. doi: 10.1002/ejhf.2073. Epub 2020 Dec 29. PMID: 33269486.
49. Januzzi JL Jr, Xu J, Li J, Shaw W, Oh R, Pfeifer M, Butler J, Sattar N, Mahaffey KW, Neal B, Hansen MK. Effects of Canagliflozin on Amino-Terminal Pro-B-Type Natriuretic Peptide: Implications for Cardiovascular Risk Reduction. *J Am Coll Cardiol*. 2020 Nov 3;76(18):2076-2085. doi: 10.1016/j.jacc.2020.09.004. PMID: 33121714.
50. Bayes-Genis A, Docherty KF, Petrie MC, Januzzi JL, Mueller C, Anderson L, Bozkurt B, Butler J, Chioncel O, Cleland JGF, Christodorescu R, Del Prato S, Gustafsson F, Lam CSP, Moura B, Pop-Busui R, Seferovic P, Volterrani M, Vaduganathan M, Metra M, Rosano G. Practical algorithms for early diagnosis of heart failure and heart stress using NT-proBNP: A clinical consensus statement from the Heart Failure Association of the ESC. *Eur J Heart Fail*. 2023 Nov;25(11):1891-1898. doi: 10.1002/ejhf.3036. Epub 2023 Sep 26. PMID: 37712339.
51. Anderson L, Bayes-Genis A, Bodegård J, Mullin K, Gustafsson S, Rosano GMC, Sundström J. Suspected de novo heart failure in outpatient care: the REVOLUTION HF study. *Eur Heart J*. 2025 Apr 22;46(16):1493-1503. doi: 10.1093/eurheartj/ehaf034. PMID: 39935142; PMCID: PMC12011520.